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EN | SEECODE Dynamic Back Massager with Heat Function

Thank you very much for purchasing our Dynamic Back Massager 2.0. With this back massager, you can enjoy a soothing massage along your spine at any time – from your lower back right up to your shoulders.

Two ergonomically shaped, rotating massage heads move gently up and down, providing a consistent and pleasantly intense massage experience. The switchable heat function further aids muscle relaxation. The device is cordless with a built-in rechargeable battery, offering maximum flexibility in everyday life.

Please read this user manual carefully before first use. Keep this user manual in a safe place. It contains important information on the safe use and care of your Dynamic Back Massager 2.0.

Intended use

The device is intended solely for the relaxing massage of the back area in a private household.

The product is intended for:

- Massage of the upper and lower back
- Massage along the spine
- Relaxation of tense muscles
- Promotion of personal well-being

The device is not intended for commercial or medical-therapeutic purposes and is not a substitute for medical treatment.

Contents

- 1 × Dynamic Back Massager 2.0
- 1 × USB-C charging cable
- 1 × User manual

The fastening strap is already attached to the device.

Technical Specifications

Product dimensions: approx. 48 × 39 × 9 cm

Weight: 1.5 kg

Operating voltage: 5V DC

Charging current: 2A

Power: 30 W

Battery: Lithium-ion battery

Battery capacity: 2,000 mAh

Charging time: approx. 4 hours

Operating time: approx. 1 hour

Automatic shut-off: after approx. 15 minutes

Heat function: approx. 42 °C

Connection: USB-C

Controls

The control panel is located on the side of the device.
All massage functions can be conveniently controlled using the four function buttons.

1 USB-C charging port

For charging the built-in battery.

2 LED indicator light

Displays the device's current charge status.

Red light: Battery is charging.

White light: Battery is fully charged.

Used for easy monitoring of the battery charge.

3 On/Off button

Press and hold the button for approx. 2 seconds to turn the device on or off.

After turning on, the massage heads automatically begin their up-and-down movement.

4 Mode button

Each press of the button changes the massage areas:

1. Full back massage
2. Upper and middle back
3. Lower and middle back

5 Speed function

The default speed is low.

Pressing the speed button once switches to high speed.

6 Heat button

The heat function is automatically activated when the device is switched on.

Press the heat button once to switch off the heat function.

Pressing it again reactivates the heat function.

Charging the battery

Fully charge the battery before using the device for the first time.

1. Connect the USB-C cable to the device's charging port.
2. Connect the cable to a 5 V / 2 A USB power supply (not included).
3. The indicator light will glow red whilst charging.
4. Once fully charged, the light will turn white.
5. Then disconnect the device from the mains.

Important:

- The device must not be used whilst charging.
- Use only the USB cable supplied.
- Using unsuitable chargers may damage the device.

Operation

Positioning the device

Place the back massager between your back and a suitable backrest, such as a chair, armchair or sofa. The integrated fastening strap can be used for additional support. Lean back comfortably and adopt a relaxed sitting position.

The intensity of the massage is determined by the amount of pressure you apply. The harder you lean against the device, the more intense the massage will feel.

Start the massage

Press and hold the On/Off button for about 2 seconds to switch the device on. Once switched on, the full-back massage starts automatically. The two massage units move gently up and down along the spine, massaging the upper, middle and lower back.

Select massage area

You can use the mode button to choose between different massage areas:

- Full back massage for a comprehensive treatment of the entire back
- Upper and middle back for a targeted massage in the shoulder and upper back area
- Lower and middle back to relax the lumbar and lower back area

This allows you to tailor the massage to your personal needs.

Adjusting the speed

The massage starts at a low speed by default. If required, you can switch to a higher massage speed by pressing the speed button to achieve a more intense treatment.

Using the heat function

The soothing heat function is automatically activated when the device is switched on. The massage heads heat up to around 42 °C, helping to loosen the muscles and promote a pleasant feeling of relaxation.

The heat function can be switched on or off at any time by pressing the heat button.

Automatic shut-off

For safety reasons, the device is equipped with an automatic shut-off function. After approximately 15 minutes of use, the device switches itself off. It can then be switched on again for further use.

Recommended usage

For an optimal massage experience, we recommend a session duration of 10 to 15 minutes. Continuous use should not exceed 30 minutes.

Regular use can contribute to a pleasant feeling of relaxation and help to relieve everyday tension in the back.

Cleaning and care

Switch off the device and disconnect it from the charger.

Clean only with a well-wrung, slightly damp cloth.

Allow to dry completely afterwards.

Do not use aggressive cleaning agents.

Never clean the device under running water.

Never immerse the device in water.

Protect from direct sunlight.

Safety instructions

WARNING

Do not use the device:

- whilst sleeping
- whilst charging
- whilst driving
- in a damp environment
- near heat sources
- if there is visible damage
- after a fall
- if moisture has entered the device

The device is intended for private use only.

Keep the device out of the reach of children.

Health advice and contraindications

This device is not suitable for:

- Pregnant women
- Children
- People with pacemakers
- People with electronic implants
- People with artificial joints
- People with varicose veins
- People with circulatory disorders
- People with reduced sensitivity to heat
- People with bruises or contusions
- People with open wounds
- People with skin irritations
- People with muscle or tissue inflammation
- People with phlebitis

If in doubt, consult your doctor before use.

Stop using the device immediately if:

- Pain occurs,
- Numbness occurs,
- Dizziness occurs,
- You feel unwell.

Disclaimer

This product is intended solely for relaxation and personal well-being.

It is not a medical device. The device is not intended for the diagnosis, treatment, cure or prevention of diseases.

The manufacturer, importer and distributor accept no liability for damage or injury resulting from improper use, failure to follow these instructions, or unauthorized modifications to the product.

Subject to technical changes and errors.

General safety instructions for non-removable batteries

- Children must only charge batteries under adult supervision. Please always keep out of reach of children.
- The electrolytes and electrolyte vapors contained in a battery are harmful to health. Avoid direct contact at all costs. If electrolytes come into contact with the skin, eyes or other parts of the body, rinse immediately with plenty of fresh water. Consult a doctor afterwards.
- Stop using or charging the battery immediately if you notice that it is swelling, leaking, becoming very hot, or emitting smoke or a strong odor.
- If you do not intend to use the device for an extended period, charge it occasionally. Connected batteries may become deeply discharged and must then no longer be used or recharged.
- Do not immerse batteries or the device in water.
- Do not short-circuit the battery (connect the positive and negative terminals).
- Do not dismantle or open batteries, as this may cause burns from the electrolytes.
- Do not expose batteries to high heat, e.g. direct sunlight, fire or similar sources. Do not place them on a radiator.
- Do not dispose of batteries by breaking them up, burning them or putting them in household waste. Please follow local regulations when disposing of them.
- Please note the legal obligation for end users to return used batteries and the option to return them free of charge to retailers.

Warranty & Guarantee

Statutory warranty

This product is covered by statutory warranty rights.

The statutory warranty period is 24 months from the date of purchase.

We provide a 6-month warranty on the battery.

You are entitled to these rights as the purchaser, regardless of any additional voluntary warranty that may be provided.

Voluntary Warranty

In addition to the statutory warranty, we provide a voluntary warranty of 6 months from the date of purchase.

This warranty is valid throughout Europe and requires the presentation of a proof of purchase.

The voluntary warranty covers material and manufacturing defects that may occur despite careful quality control prior to dispatch.

The following are specifically excluded from the voluntary warranty:

- Alterations or modifications to the product without our consent
- Use of the product in a manner that does not correspond to its intended use

Procedure in the event of a warranty or guarantee claim

Please return the complete product in secure packaging to the point of sale where it was purchased.

Please enclose the following with your return:

- a description of the defect found,
- the original proof of purchase as evidence for the warranty or guarantee.

The seller will examine the reported defect and reserves the right to rectify the fault before issuing any refund of the purchase price.

In the event of a valid warranty or guarantee claim, all necessary transport costs will be reimbursed.

The exercise of statutory rights in respect of defects is free of charge for the buyer.

General information

© Copyright

These instructions are protected by copyright. Reproduction or illustration, even in modified form, is permitted only with the written consent of the manufacturer.

CE Declaration of Conformity

Mobiset GmbH hereby declares that this product complies with the essential requirements of the applicable European directives.

Environmental Protection Information

Date of placing on the market: 01/10/2026

Disposal of the product

Please dispose of the product in an environmentally responsible manner in accordance with applicable local regulations:

- This product is an electrical appliance containing a built-in lithium-ion battery and must not be disposed of with household waste.
- Take the appliance to an official collection point for waste electrical and electronic equipment or a recycling center.
- Rechargeable and non-rechargeable batteries must not be disposed of with household waste, as they may contain harmful substances. If possible, remove the battery correctly or hand in the device in its entirety at a suitable collection point.
- Electrical and electronic devices contain valuable raw materials that can be reused through recycling.

Disposal of packaging

Packaging materials (e.g. cardboard and plastic film) should be separated and disposed of in accordance with local recycling guidelines.

By disposing of the product and packaging correctly, you are making an important contribution to environmental protection and the sustainable use of resources.

Please contact your local waste management authority for information on the applicable disposal and recycling options.

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